DeVilbiss® iFill™ Oxygen Cylinder

DeVilbiss® iFill™ Oxygen Cylinder

PD1000A Series iFill Oxygen Cylinder w/Integrated PulseDose® Conserving Device

535-CF Series iFill Oxygen Cylinder w/Continuous Flow Regulator

CAUTION-Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

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Improving People's Lives

TABLE OF CONTENTS

Important Safeguards	4
Physician's Notes	5
Introduction	5
Use Times	6
Important Parts of Your iFill Oxygen Cylinder.....	7
Cylinder Labeling Verbiage.....	8
Operating Instructions	9
Inserting Batteries into the iFill Oxygen Cylinder w/PulseDose	9
Using Your Oxygen iFill Cylinder w/PulseDose	9
Using Your Oxygen iFill Cylinder w/Continuous Flow Regulator	9
Typical Questions and Answers about PulseDose	10
Care & Maintenance	11
Troubleshooting.....	11
Specifications.....	12
iFill Oxygen Cylinder w/Integrated PulseDose	12
iFill Oxygen Cylinder w/Continuous Flow Regulator	12
All iFill Oxygen Cylinders	12

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IMPORTANT SAFEGUARDS

The information contained in this guide is intended to assist in the safe operation of the equipment and to ensure maximum benefit is achieved.

This product should only be filled by the DeVilbiss iFill Personal Oxygen Station.

READ ALL INSTRUCTIONS BEFORE USING. SAVE THESE INSTRUCTIONS.

Oxygen supplied by this equipment is not to be considered life-supporting and must not supply anything other than 93% ± 3% oxygen.

When using electrical products, basic safety precautions should always be followed. Read all instructions before using.

Important information is highlighted by these terms:

DANGER– Urgent safety information for hazards that will cause serious injury or death.

WARNING– Important safety information for hazards that might cause serious injury.

CAUTION– Information for preventing damage to the product.

NOTE– Information to which you should pay special attention.

 **ATTENTION**–Consult accompanying documents.

 **DANGER– NO SMOKING**

 Type BF applied part

DANGER

To reduce the risk of fire, burns, or injury to persons:

Oxygen, though non-flammable, vigorously supports and accelerates burning of any flammable material. If you know or suspect oxygen has escaped other than through normal operation, open doors and windows to ventilate the area.

1. **DO NOT SMOKE WHILE USING YOUR DEVILBISS OXYGEN EQUIPMENT.** Keep matches, cigarettes, burning tobacco, or candles away from the area where the system is being stored or operated.
2. Avoid creation of any spark near oxygen equipment. This includes sparks from static electricity created by any type of friction.
3. Keep the equipment at least seven feet away from radios, television sets, window air conditioners, fans, electric razors, hair dryers, and all other electrical appliances.
4. Keep the equipment at least five feet away from heat sources, electric or gas heaters of any kind, fireplaces, or stoves.
5. Keep all flammable materials or petroleum-based products away from the equipment.
6. Never attempt to lubricate the equipment.
7. Never use aerosol sprays near the equipment.

To prevent high concentrations of oxygen:

1. Keep the equipment in a well-ventilated area.
2. Do not carry equipment under a coat or any form of clothing.
3. Turn off oxygen supply by turning the rotary selector to the Off position when not in use.

WARNING

To reduce the risk of injury:

1. Keep all units away from children. Do not allow unauthorized or untrained individuals to operate the equipment. Never tamper with or try to repair the equipment yourself. If you have any questions or suspect your equipment is not operating properly, contact your oxygen provider.
2. Do not immerse in liquids or subject device to harsh conditions.
3. Do not use with other equipment (i.e. humidifier, nebulizer, etc.) when in PulseDose delivery mode.

PHYSICIAN'S NOTES

1. Use only continuous flow iFill oxygen cylinders (535-CF series) with patients who breathe below 6 Breaths Per Minute (BPM) or above 40 BPM.
2. Use only continuous flow iFill oxygen cylinders (535-CF series) with patients who consistently fail to trigger equipment (i.e. mouth breathing with closed soft palates).
3. Verify patient is getting adequate PaO₂ or SaO₂ levels in PulseDose delivery.
4. Use only standard nasal cannula with PulseDose delivery. Do not use pediatric (low-flow) nasal cannula or mask with PulseDose delivery.
5. A mask or any nasal cannula can be used with continuous flow delivery.
6. The PD1000A series electronic conserving device contains a preset 2 LPM continuous flow backup. The service manual describes how to change the cannula fitting to obtain a 3, 4, 5 or 6 LPM continuous flow backup.

INTRODUCTION

INDICATIONS FOR USE (PD1000A SERIES)

The DeVilbiss iFill Oxygen Cylinders w/integrated PD1000A are intended as a delivery device for 93% ± 3% oxygen from high-pressure oxygen cylinders. This is an ambulatory device, which allows patients to ambulate longer than they would with a continuous flow regulator on the same cylinder.

HOW PULSED DOSE WORKS

PulseDose dramatically extends the use time from a supply of oxygen, offering increased mobility with improved comfort and increased efficiency. The reliability and safety of PulseDose oxygen delivery has been proven effective in clinical testing as well as through independent tests performed by physicians and respiratory therapists.

What is PulseDose? The concept is based on the fact that the normal breathing pattern is inhalation for 1/3 of the time, and exhalation about 2/3 of the time. At 20 BPM the oxygen delivered in continuous flow, assuming inspiration is 1/3 of the breathing cycle, would be 16.5 cc/LPM. As a result, PulseDose extends the use time of an oxygen system by an average of 3:1. PulseDose senses the start of inhalation and instantly releases a short “pulsed” dose at the very beginning of the breathing cycle. Since all of the “pulsed” oxygen finds its way deep into the lungs, less oxygen is required to accomplish the same effect than with traditional continuous flow oxygen systems. This means that a PulseDose oxygen system will last two to four times longer than a continuous flow oxygen system, yet still provide the same therapeutic benefit.

Because oxygen is released only during inhalation, the constant flow of oxygen into the nostrils is eliminated. Many users find PulseDose oxygen delivery more comfortable than continuous flow delivery systems. The short “pulse” of oxygen delivered during inhalation is almost undetectable, and the humidity in the room air helps maintain a normal level of moisture in the nasal cavity. This greatly reduces the discomfort of dehydration associated with a continuous flow oxygen system.

USE TIMES

Because PulseDose responds to each individual's breathing patterns, the use time will vary for each individual depending on the PulseDose prescription rate and the breath rate. The following chart shows the theoretical ambulatory ranges for DeVilbiss PulseDose products.

NOTE—All ambulatory ranges are calculated assuming a breath rate of 20 breaths per minute in PulseDose (PD) mode.

USE TIMES SHOWN IN HOURS									
Delivered Volume cc's:	16.5	24.75	33	41.25	49.5	66	82.5	99	
Flow Rate:	1	1.5	2	2.5	3	4	5	6	Mode
iFill M-4 Cylinder 113 Gaseous Liters	1.9	1.3	.9	.7	.6	.5	.4	.3	CF
	5.7	3.8	2.9	2.3	1.9	1.4	1.1	.9	PD
iFill M-6 Cylinder 164 Gaseous Liters	2.7	1.8	1.4	1.1	.9	.7	.6	.4	CF
	8.3	5.5	4.1	3.3	2.8	2.1	1.7	1.4	PD
iFill ML-6 Cylinder 170 Gaseous Liters	2.8	1.9	1.4	1.1	.9	.7	.6	.5	CF
	8.6	5.7	4.3	3.4	2.9	2.1	1.7	1.4	PD
iFill C Cylinder 240 Gaseous Liters	4.0	2.7	2.0	1.6	1.3	1.0	.8	.7	CF
	12.1	8.1	6.1	4.9	4.0	3.0	2.4	2.0	PD
iFill D Cylinder 415 Gaseous Liters	6.9	4.6	3.5	2.8	2.3	1.7	1.4	1.2	CF
	21.0	14.0	10.5	8.4	7.0	5.2	4.2	3.5	PD
iFill E Cylinder 682 Gaseous Liters	11.4	7.6	5.7	4.6	3.8	2.8	2.3	1.9	CF
	34.4	23.0	17.2	13.8	11.5	8.6	6.9	5.8	PD


CF=Continuous Flow

PD=PulseDose

Specifications subject to change without notice. This chart is intended to be used only as a guide.

Cylinders vary in gaseous liter capacity by manufacturer which may result in varying use times.

IMPORTANT PARTS OF YOUR IFILL OXYGEN CYLINDER

PD1000A SERIES

1. Cannula Fitting– Use this fitting to attach the cannula to your iFill cylinder.
2. Rotary Selector– Use this to select your prescription setting.
3. Oxygen Fill Connection – Use this connection to attach the cylinder to the DeVilbiss iFill Personal Oxygen Station. **NOTE**–Your iFill oxygen cylinder can only be filled with the DeVilbiss iFill Personal Oxygen Station.
4. Oxygen Contents Gauge– Indicates the remaining pressure in the oxygen cylinder. When this gauge falls into the red section, you should switch to a new cylinder.
5. Hydrostatic Test Date - Contact your homecare provider or the cylinder manufacturer for details.
6. High-Pressure Oxygen Cylinder – Portable oxygen delivery tank with labeling. Refer to next page for label verbiage.
7. • PulseDose Indicator– Either green or red light illuminates each time the unit pulses oxygen.
• Normal Battery Indicator– A flashing green light indicates that there is sufficient battery power.
8. • Low Battery Indicator– A flashing red light indicates that there are 4 - 8 hours of battery life remaining. The low battery indication time may be reduced when using NiMH batteries.

- Change Battery Indicator– A constant red light indicates that the battery should be changed immediately. Unit can be used only in continuous flow mode until new batteries are installed.
9. Battery Door– Use only a standard “AA” Alkaline or NiMH batteries.

535-CF SERIES

1. Cannula Fitting– Use this fitting to attach the cannula to your iFill cylinder.
2. Rotary Selector– Use this to select your prescription setting.
3. Oxygen Fill Connection – Use this connection to attach the cylinder to the DeVilbiss iFill Personal Oxygen Station. **NOTE**–Your iFill oxygen cylinder can only be filled with the DeVilbiss iFill Personal Oxygen Station.
4. Oxygen Contents Gauge– Indicates the remaining pressure in the oxygen cylinder. When this gauge falls into the red section, you should switch to a new cylinder.
5. Hydrostatic Test Date - Contact your homecare provider or the cylinder manufacturer for details.
6. High-Pressure Oxygen Cylinder – Portable oxygen delivery tank with labeling. Refer to next page for label verbiage.

CYLINDER LABELING VERBIAGE

THIS CYLINDER IS TO BE USED WITH THE DeVILBISS iFILL PERSONAL OXYGEN STATION. THIS CYLINDER CONTAINS 93+/-3% OXYGEN PRODUCED BY THE DeVILBISS iFILL PERSONAL OXYGEN STATION USING A PRESSURE SWING ADSORPTION PROCESS.

OXYGEN CYLINDERS FILLED BY THE DeVILBISS iFILL PERSONAL OXYGEN STATION SHALL BE USED FOR PERSONAL USE ONLY. "NOT TO BE FILLED FOR RESALE"

WARNING: HIGH PRESSURE OXIDIZING GAS VIGOROUSLY ACCELERATES COMBUSTION.

CAUTION: For medical applications. Use only as directed by a licensed practitioner. Uninterrupted use of high concentrations of oxygen over a long duration, without monitoring its effect on oxygen content of arterial blood, may be harmful.

CAUTION: DO NOT HANDLE CYLINDER OR USE CONTENTS UNTIL YOU ARE TRAINED TO USE CYLINDER AND CONTENTS, INCLUDING EMERGENCY PROCEDURES.

CAUTION: IMPROPER FILLING OR ABUSE OF THIS CYLINDER, OR FAILURE TO HEED THIS CAUTION MAY CAUSE SERIOUS INJURY OR LOSS OF LIFE. DO NOT ALTER OR MODIFY CYLINDER OR RELATED COMPONENTS. KEEP OUT OF THE REACH OF CHILDREN. Secure cylinder during storage and use. No smoking in cylinder area. Turn the regulator or conserver off after each use and when empty. Use in accordance with the DeVilbiss iFILL Personal Oxygen Station instruction guide. Cylinder temperature should not exceed 130°F (54°C). Keep away from heat, flame and spark. Do not drop.

"Use No Oil or Grease". KEEP CYLINDER, ALL EQUIPMENT AND CONNECTIONS FREE OF OIL AND GREASE TO AVOID VIOLENT IGNITION.

TRANSFILLING OF THIS GAS IS PERFORMED BY THE DeVILBISS iFILL PERSONAL OXYGEN STATION ONLY.

Cylinders that have been refinished and/or exposed to elevated temperatures must be hydrostatically tested before filling. Do not use caustic paint stripper or corrosive cleaners. Aluminum cylinders subject to the action of fire or heated to temperatures in excess of 350°F (177°C) must be withdrawn from service and condemned by trained personnel. Valve and safety relief devices must be removed and replaced by trained and authorized personnel. **Do not alter or modify this cylinder in any way.**

DO NOT REMOVE THIS PRODUCT LABEL
Not for Commercial Resale

Sunrise Medical
Respiratory Division
100 DeVilbiss Drive
Somerset PA 15501

**WARNING****ALUMINUM HIGH PRESSURE GAS CYLINDER**

Improper use, filling, storage, or disposal may cause personal injury, death, or property damage. Do not alter or modify this cylinder or valve in any way. Do not use caustic or corrosive cleaners. Never fill this cylinder unless it has been hydrostatically tested. Only trained personnel should fill or service the cylinder. Only trained personnel should replace valves and pressure relief devices and only with the complete assemblies containing the proper burst disc rating supplied by the valve manufacturer. Always secure cylinder in a cool dry area, out of reach of children. Do not expose filled cylinder to temperatures above 130°F (54°C). Cylinders exposed to fire or heat in excess of 350°F (177°C) must be condemned. Cylinders refinished or subjected to elevated temperatures must be hydrostatically tested prior to refilling. Do not remove, alter, or obscure this warning label. Proper cleaning of cylinder, valve, and components is required before filling with oxidizing gasses such as oxygen and nitrous oxide.

Service, maintain and inspect cylinders in accordance with C.G.A. pamphlet C-1, C-6.1, G-4, G-4.1, G-4.3, P-1, P-2, P-2.5, P-14, and G-8.2 available from the CGA at ph. 703-788-2700 or web: www.cganet.com

MADE IN U.S.A.

PIPER METAL FORMING CORPORATION

NEW ALBANY, MS

www.piperimpact.com

OPERATING INSTRUCTIONS

INSERTING BATTERIES INTO THE IFILL OXYGEN CYLINDER W/PULSEDOSE

1. Open the battery door.
2. Insert 2 "AA" Alkaline or NiMH batteries.
3. Close the battery door.

NOTE— When changing batteries, first turn the rotary selector to the "OFF" position.

USING YOUR IFILL OXYGEN CYLINDER W/PULSEDOSE

1. Attach the standard nasal cannula to the cannula fitting and to your nose and face. Oxygen tubing up to 35 feet long may be used in PulseDose delivery mode.

⚠ WARNING

To prevent injury from cylinders tipping over, do not use cannula tubing lengths over 10 feet with small compressed oxygen cylinders. Unattended cylinders should be secured in a cylinder stand.

2. Turn the rotary selector to the prescribed flow setting.

NOTE— When this rotary switch is set to "OFF," the unit is not using battery power and will not pulse. When the selector is set to one of the numbers, the unit is on and awaiting inspiration through the nasal cannula at which time it will dose on every breath. The volume of the oxygen delivered varies according to which prescription flow setting is chosen. The final setting on the rotary switch is "CF"; this is the continuous flow position. In this position oxygen will flow from the cannula fitting at the default continuous flow rate.

NOTE— Continuous flow mode is not powered by the batteries and can be used regardless of the battery level. In the event of a device failure or dead batteries, the user must manually switch the unit to continuous flow for delivery of oxygen. The device will not automatically switch to continuous flow. The oxygen cylinder will not last as long in continuous flow mode as it would in PulseDose mode. Unless there is a problem with the unit, such as dead batteries, the unit should be used in PulseDose mode.

3. Breathe normally, the conserving device will deliver a bolus of oxygen at the leading edge of inspiration on every breath up to 40 breaths per minute.

4. When you are finished using the iFill oxygen cylinder, turn the rotary selector to the "OFF" position.

NOTE— A mask should not be used in the PulseDose delivery mode as it may not fit to the face well enough to allow the conserving device to sense inhalation efforts. Also, the therapeutic effect of PulseDose would not be realized, as the dose of oxygen would be diluted in the mask prior to inhalation.

NOTE— A pediatric or low-flow cannula should not be used in PulseDose delivery mode. The reduced diameter of the cannula causes too much back pressure and will affect the oxygen volume delivered.

NOTE— PulseDose delivers oxygen in a very short "puff." It does not deliver oxygen during the entire inhalation. The length of time that PulseDose delivers oxygen will not vary from breath to breath. The time is set in correlation to the oxygen dosage set on the conserving device (patient's prescription setting).

NOTE— PulseDose is designed to prevent the delivery of pulses more than every 1 1/2 seconds. If the breath rate is greater than 40 BPM, this feature prevents delivery of excessive oxygen by not dosing on every breath.

NOTE— If using NiMH batteries, carefully monitor when the low battery indicator flashes red. It is recommended that a spare fully-charged set of alkaline batteries be kept in reserve and installed when the change battery indicator is constant red.

ALERTS (PD1000A-CAN-M4, -M6, -ML6, -C, -D, AND -E)



No Inspiration Alert – If the unit is on and inspiration has not been sensed for 30 seconds (after at least one breath has been sensed), the audible alert will activate every three seconds for one minute and then auto shut-off.

USING YOUR IFILL OXYGEN CYLINDER W/CONTINUOUS FLOW REGULATOR

1. Attach the nasal cannula (or mask) to the cannula fitting and to your nose and face.
2. Turn the rotary selector to the prescribed flow setting. When the selector is set to one of the numbers, the unit is on.
3. Breathe normally.
4. When you are finished using the iFill oxygen cylinder, turn the rotary selector to the "OFF" position.

⚠ WARNING

To prevent injury from cylinders tipping over, do not use cannula tubing lengths over 10 feet with small compressed oxygen cylinders. Unattended cylinders should be secured in a cylinder stand.

TYPICAL QUESTIONS AND ANSWERS ABOUT PULSEDOSE**Q. How does PulseDose work? How does it know when I'm inhaling?**

A. When inhaling, your diaphragm moves down and causes a drop in pressure in the lungs. Air flows in through the nose and mouth to equalize the pressure. This negative pressure is also present at the nose and mouth during inhalation. This pressure signal travels through the nasal cannula to a pressure sensor in the PulseDose conserving device. An electronic circuit then opens an electrical valve to deliver a precisely metered dose of oxygen. When the valve is closed, the sensor is ready to detect the next inhalation.

Q. The pulse seems so short. Am I really getting enough oxygen?

A. Yes. PulseDose delivers a precise burst of oxygen at a relatively high flow rate at the leading edge of each inhalation. This assures that the oxygen delivered flows deep into the lungs for maximum benefit. PulseDose requires less oxygen to deliver the same therapeutic benefit as continuous flow oxygen delivery.

Q. I can't hear the pulse. Is PulseDose working?

A. If the pulse can't be heard, simply look at the green PulseDose indicator to see that the device is being triggered by inhalation. For further assurance, hold the end of the cannula in front of your lips while inhaling through your mouth and feel the pulse. PulseDose does not monitor the supply of oxygen. Remember to check the oxygen contents gauge periodically to verify that there is an adequate oxygen supply. If the oxygen supply runs out, the green PulseDose indicator light will continue to illuminate, indicating that the conserving device is being triggered by inhalation.

Q. Why can't I use a cannula which is longer than 35 feet with PulseDose?

A. The PulseDose triggering is not significantly affected by the cannula length, but the delivery of oxygen is affected. If the cannula is longer than 35 feet, the pulse of oxygen is delayed. Remember the therapeutic moment during the inhalation cycle. If the oxygen is not delivered during this time, the benefits will not be realized.

Q. I've always used humidifiers with oxygen. Should I use a humidifier with PulseDose?

A. No. PulseDose is not able to sense inhalation through the water in the humidifier. Also, many patients find that humidification is not necessary with PulseDose. They find that PulseDose improves comfort because it delivers a very small amount of oxygen during the early part of inhalation, while the rest of the inhalation is composed of normal room air.

Q. When I'm breathing faster, I don't get a pulse with each breath. Don't I need a dose every time?

A. Because PulseDose breathes with the patient, it has an upper limit (40 Breaths Per Minute) that keeps you from getting too much oxygen. When breathing slowly, you receive a dose with every breath. As breath rate increases (up to 40 BPM) PulseDose still delivers a dose with every breath. At this point, you are getting more oxygen per minute because each pulse delivers the same amount of oxygen with each breath while the number of breaths has increased. With continuous flow oxygen, the oxygen delivered is constant. As you breathe faster, the enrichment of inhalations actually decreases because each breath is being diluted with a greater amount of room air.

CARE & MAINTENANCE

The DeVilbiss iFill Oxygen Cylinders should be kept clean and free from moisture and dust. The cylinders should be protected from extreme temperatures. Clean the cylinders periodically by wiping with a dry, lint-free cloth. Do not clean with a solvent based cleaning solution. Avoid dropping the cylinder or placing it in a position where it could topple or fall since this can damage the device. Whenever possible, use a padded carrying bag such as the Hideaway or Walkabout bags to carry the cylinder. This will help to protect it in case of a fall. Avoid getting fluids or debris such as sand or dirt inside the oxygen cylinder and/or oxygen connector.

STORAGE AND HANDLING WARNINGS

Do not place oxygen cylinders in unventilated spaces such as car trunks. Excessive heat can make the relief valve suddenly and quickly discharge the cylinder contents, possibly making it a projectile and greatly increasing the oxygen level in unventilated spaces.

Do not leave oxygen cylinders in the cabs of vehicles without ventilation. If a cylinder leaks, a spark could start a fire causing serious injury or death.

Remove cylinders from the vehicle when the destination is reached.

Be sure to secure cylinders from movement during transport.

TROUBLESHOOTING

WARNING

Do not attempt to open the device for maintenance or repair. The device contains no user-serviceable parts. Contact your oxygen provider if service is required.

Symptoms	Possible Causes	Remedies
Oxygen is not being delivered even though the PulseDose indicator is flashing every time I breathe.	1. Oxygen supply is empty.	1. Check contents indicator on the device. If empty, switch cylinders.
Use times are different from those stated in the literature.	1. PulseDose responds to your breath rate. Your breath rate may vary, which causes the operation time to vary. 2. Leak in system	1. PulseDose is operating correctly. 2. Contact your Sunrise Medical provider.
PulseDose will not pulse	1. Cannula is not adjusted properly. 2. Unit is not turned on. 3. Batteries discharged or not installed. 4. Mouth breathing with closed soft palates. 5. Unit did not reset while changing batteries (red light stays on).	1. Check all cannula connections to make sure they are tight, and adjust the cannula to fit comfortably in your nose. Ensure tubing is not kinked. 2. Turn the rotary selector to the appropriate setting. 3. Install new batteries. 4. Breathe through the nose (cannula). 5. Turn the unit off and back on using the rotary selector.
PulseDose works fine for a couple of minutes, then sensitivity seems to drift and may stop working altogether.	1. Using pediatric cannula or any cannula that restricts continuous flow capacity of 10 lpm.	1. Replace with standard nasal cannula.

SPECIFICATIONS

IFILL OXYGEN CYLINDER W/INTEGRATED PULSED DOSE

iFill Oxygen Cylinder w/Integrated PulseDose Regulator (US)		iFill Oxygen Cylinder w/Integrated PulseDose Regulator (CAN)	
Weight			
PD1000A-M4	3.5 lbs (1.59 kg)	PD1000A-CAN-M4	3.5 lbs (1.59 kg)
PD1000A-M6	4.1 lbs (1.86 kg)	PD1000A-CAN-M6	4.1 lbs (1.86 kg)
PD1000A-ML6	4.8 lbs (2.18 kg)	PD1000A-CAN-ML6	4.8 lbs (2.18 kg)
PD1000A-C	5.6 lbs (2.54 kg)	PD1000A-CAN-C	5.6 lbs (2.54 kg)
PD1000A-D	7.2 lbs (3.27 kg)	PD1000A-CAN-D	7.2 lbs (3.27 kg)
PD1000A-E	9.8 lbs (4.45 kg)	PD1000A-CAN-E	9.8 lbs (4.45 kg)
Dimensions			
PD1000A-M4	3.2"D x 12.4"L (81mm x 315mm)	PD1000A-CAN-M4	3.2"D x 12.4"L (81mm x 315mm)
PD1000A-M6	3.2"D x 15.59"L (81mm x 396mm)	PD1000A-CAN-M6	3.2"D x 15.59"L (81mm x 396mm)
PD1000A-ML6	4.38"D x 11.68"L (111mm x 297mm)	PD1000A-CAN-ML6	4.38"D x 11.68"L (111mm x 297mm)
PD1000A-C	4.38"D x 14.88"L (111mm x 378mm)	PD1000A-CAN-C	4.38"D x 14.88"L (111mm x 378mm)
PD1000A-D	4.38"D x 20.51"L (111mm x 521mm)	PD1000A-CAN-D	4.38"D x 20.51"L (111mm x 521mm)
PD1000A-E	4.38"D x 29.63"L (111mm x 753mm)	PD1000A-CAN-E	4.38"D x 29.63"L (111mm x 753mm)
Power Supply(2) Standard "AA" alkaline or NiMH batteries.		
Power RequirementsBatteries other than alkaline or NiMH are not recommended due to the capacity needed for operation and battery life of the unit.		
Degree of Protection Against Electric ShockTYPE BF applied part		
Modes of OperationContinuous / Pulsed		
Approval Body And Safety StandardIEC 601-1; CAN/CSA-C22.2 No 601.1-M90		
US Patents4,519,387; 5,755,224; 4,457,303		

IFILL OXYGEN CYLINDER W/CONTINUOUS FLOW REGULATOR

iFill Oxygen Cylinder w/Continuous Flow Regulator (US)		iFill Oxygen Cylinder w/Continuous Flow Regulator (CAN)	
Weight			
535D-M6-CF	3.9 lbs (1.77 kg)	535C-M6-CF	3.9 lbs (1.77 kg)
535D-ML6-CF	4.6 lbs (2.09 kg)	535C-ML6-CF	4.6 lbs (2.09 kg)
535D-C-CF	5.4 lbs (2.45 kg)	535C-C-CF	5.4 lbs (2.45 kg)
535D-D-CF	7.0 lbs (3.18 kg)	535C-D-CF	7.0 lbs (3.18 kg)
535D-E-CF	9.6 lbs (4.35 kg)	535C-E-CF	9.6 lbs (4.35 kg)
Dimensions			
535D-M6-CF	3.2"D x 14.9"L (81mm x 378mm)	535C-M6-CF	3.2"D x 14.9"L (81mm x 378mm)
535D-ML6-CF	4.38"D x 11.0"L (111mm x 279mm)	535C-ML6-CF	4.38"D x 11.0"L (111mm x 279mm)
535D-C-CF	4.38"D x 14.2"L (111mm x 361mm)	535C-C-CF	4.38"D x 14.2"L (111mm x 361mm)
535D-D-CF	4.38"D x 19.82"L (111mm x 503mm)	535C-D-CF	4.38"D x 19.82"L (111mm x 503mm)
535D-E-CF	4.38"D x 28.94"L (111mm x 735mm)	535C-E-CF	4.38"D x 28.94"L (111mm x 735mm)

ALL IFILL OXYGEN CYLINDERS

Operating Temperature Range 41° to 104°F (5° to 40°C)
Operating Pressure Range 500 to 2250 PSIG (34 to 155 bar) tank pressure
Operating Atmospheric Conditions 500 to 1020 millibar
Operating Humidity Range 0 to 95% R.H., non-condensing
Storage and Transportation Temperature Range -4° to 130°F (-20° to 54°C)
Storage and Transportation Humidity Range Up to 95% R.H., non-condensing
Degree of Protection Against Ingress of Liquids NONE
Safety Standard meets DOT-3AL

Specifications subject to change without notice.

NOTE- Degradation of performance may occur if unit is operated outside of specified operating parameters.

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Sunrise Medical Pty. Limited